# PikoButton

# SE PikoButton

## Beskrivning

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# GB PikoButton

## **Description**

PikoButtons are control switches, designed to be very durable but still easy and pleasant to use. Because of PikoButtons' sensitivity they are suitable for controlling the environment control systems, communication devices, assistant alarms, etc.

PikoButtons are available in several models and colours which allows them to be used in several different combinations.

The activation force of the buttons is very small, only 75 grams for FingerButton, FingerDual and PikoButton30. 125 grams goes for the bigger PikoButton50. The same activation force can be used at the edge of the switch.

The mounting of the buttons is very easy: PikoButton has one M4-hole in the middle of the bottom. Finger-buttons has Velcro to be adjusted for use over two fingers.

## **Switch care**

- Wrap the cable separately
- Use a slightly damp cloth with water or alcohol
- Contact your supplier if you have problems with the switch

### **Do not:**

- wrap the cable around the switch
- pull on the cable of the switch
- immerse the switch in water
- use solvents to clean your switch
- overload your switch connection
- expose switch to extreme heat or cold

## **Technical information**

Note: The guarantee expires when the switch is opened. The manufacturer is not responsible for the damage that might occur in the opening process.

For a summary of the various models and technical information, see last page.

## **Recycling Waste Electrical and Electronic Equipment**

PikoButton contains electronic components that must be recycled. When the product is disposed of, hand it in to a recycling centre for waste electrical of electrical and electronic equipment.

## PikoButton Technical information

### PikoButton30

Diameter: 30 mm  
Höjd/Height: 15 mm  
Manöverkraft/Operating force: 75 gr  
Vikt/Weight: 40 gr  
Max belastning/Max load: 12V/50 mA  
Max/min temp: -30 – +80°C  
Kabellängd/Cable length: 1,5 m  
Kontakt/Plug: 3,5 mm

Art:	Färg/Colour
451000	Röd/Red
451001	Grön/Green
451002	Blå/Blue
451003	Gul/Yellow
421004	Orange
451005	Svart/Black
451006	Vit/White
451007	Guld/Gold
451008	Silver
451020	Röd vattentålig/ Red water resistant
451021	Blå vattentålig/ Blue water resistant

### PikoButton50

Diameter: 50 mm  
Höjd/Height: 15 mm  
Manöverkraft/Operating force: 125 gr  
Vikt/Weight: 60 gr  
Max belastning/Max load: 12/50 mA  
Max/min temp: -30 – +80°C  
Kabellängd/Cable length: 1,5 m  
Kontakt/Plug: 3,5 mm

Art:	Färg/Colour
451500	Röd/Red
451501	Grön/Green
451502	Blå/Blue
451503	Gul/Yellow
451504	Orange
451505	Svart/Black
451506	Vit/White
451507	Guld/Gold
451508	Silver
451524	Blå vattentålig/ Blue water resistant
451525	Röd vattentålig/ Red water resistant

### 451510 FingerButton

Manöverkraft/Operating force: 75 gr  
Max belast/Max load: 12V/50 mA  
Max/min temp: -30 – +80°C  
Kabellängd/Cable length: 1,5 m  
Kontakt/Plug: 3,5 mm

### 451511 FingerDual

Manöverkraft/Operating force: 75 gr  
Max belast/Max load: 12V/50 mA  
Max/min temp: -30 – +80°C  
Kabellängd/Cable length: 1,5 m  
Kontakt/Plug: 3,5 mm stereo