

RELEASE NOTES FOR ROLLTALK

Release notes for Rolltalk Designer 6.7.2 release packages

769: +46 required in Sweden

Swedish country code no longer required, region code working

2331: Screen keyboard window is maximized

Fixed problem with maximizing customized window when fixed windows size is defined.

3095: Save for Micro Rolltalk should be visible in Simple Menu

Added Micro Rolltalk menus to the simple menu level.

3096: Text placement should be set to above picture by default

Changed default text placement on fields to over picture.

3112: COM port search on Light V1 stops on COM1

Fixed 'freeze' when searching for comport on some computers.

3019: Default day colors in Sweden and Denmark

Fixed default day colors in Danish version.

3138: Dwell animation not visible in entire field

Fixed update problems with square dwell animation.

3167: Runtime error with wireless switches and activity light

Fixed stability problem related to wireless switches and activity light.

3241: Capital letters and Capslock goes out of sync

Fixed out of synch problem with Capital letters on screen keyboards when page was defined to show only Capitals.

2881: Gewa infrared signals updated

Added new page library for Gewa Ir phone.

2969: Text-to-speech engine names refreshed

Added automatic update of speech engine names that changed with new Acapela voices in 6.7.

3092: Field frame is imported to paint editor

Removed frame from picture when field content is edited by paint.

3105: COM ports move to different locations

Fixed stability problem with USB ports in Windows Vista and Windows 7, that changes during runtime.

3062: Error in Micro Rolltalk keyboards

Fixed problem with macros on pages defined to have Capital letters

3329: Implement Rasmus for Micro Rolltalk

Rasmus is now a part of the 6.7.2 release cycle.

3330: Babiledll must be updated to support Acapela 7.7

Micro Rolltalk now supports Rasmus, and Acapela has fixed some Norwegian issues