

Information om vilka bild, ljud och videoformat som Programsnickaren stödjer

	PC	Micro
Ljud	WAV WMA MPA MP2 OGG MP3 AU AIF IFF SND MID IDI RMI	WAV
Bild	Raster Color and Grayscale ABIC - ABIC Format AFP - AFP Format ANI - Windows Animated Cursor AVI - Windows AVI Format BMP - Windows Bitmap Formats CIN - Cineon Format CLP - Microsoft Windows Clipboard Format CMP - LEAD Compressed CMW - LEAD Wavelet Compressed CRW - Canon RAW Format CUR - Windows Cursors CUT - Dr. Halo Format DCR - Kodak Professional Digital Camera Format DCS - Kodak Professional Digital Camera System Format DCX - DCX Format DICOM - DICOM Format DJVU - DjVu Format ECW - Enhanced Compressed Wavelet Format EMF - Windows Metafile Formats EPS - Encapsulated PostScript Format EXIF - Exif Formats FIT - FITS Format FLC - Flic Animation Format FPX - Kodak Format GIF - Compuserve GIF Format ICO - Windows Icons IFF - Interchange File Formats JBG - JBIG Format	Png vid överföring till Micro. Micro stödjer även jpg-filer som är tagna med kameran i Micro.

JPG - JPEG Compressed
JP2 - JPEG 2000 Compressed
KDC - Kodak Digital Camera Format
MPG - MPEG Video Format
MRC - Mixed Raster Content
MRC - LEADTOOLS Mixed Raster Content
MRC - LEADTOOLS PDF with Mixed Raster Content
PBM - Portable Bitmap Utilities Format
PCD - Kodak Format
PCT - Macintosh Pict Format
PCX - PCX Format
PDF - Portable Document Format
PNG - Portable Network Graphics Format
PS - PostScript Document Format
PSD - PhotoShop 3.0 Format
PSP - Paint Shop Pro Format
RAS - SUN Raster Format
RTF - Rich Text Format
SCT - Scitex Continuous Tone Format
SFF - Structured Fax File Format
SGI - Silicon Graphics Image Format
SID - MrSID (32-bit) Format
TGA - Truevision TARGA Format
TIFF - Tagged Image File Format
TIFX - Xerox Internet Fax File Format
TXT - Text Format
WBMP - Wireless Bitmap Format
WMF - Windows Metafile Formats
WPG - WordPerfect Format
XPM - XPicMap
XWD - X WindowDump

Raster Bitonal (1-Bit)

ABC - ABC Format
AWD - Microsoft FAX Format
CAL - CALS Raster Format
CMP - LEAD 1-Bit Format
FAX - Raw Fax Format
FAX - LaserView LaserData Format
ICA - Image Object Content Architecture
(IOCA/MODCA)
IMG - GEM Image Format
ITG - Intergraph Format
JBIG2 - JBIG2 Format
MAC - MacPaint Format
MSP - Microsoft Paint
SMP - Xionics Format
TIFF - Tagged Image File Format / CCITT
WFX - WinFax Format
XBM - XBitMap Format

Vector

CGM - Computer Graphics Metafile



ABILIA

GEWA AB | FALCK IGEL

	<p> CMX - Corel Presentation Exchange Format DGN - Intergraph Format DRW - Micrografx Designer DRaWing DWF - DWF Format DWG - DWG Format DXF - Drawing Interchange Format EMF - Windows Metafile Formats GBR - Gerber Format NAP - NAPLPS Format (NAP) NITF - NITF Format PCL - Printer Command Language Format PCT - Macintosh Pict Format PDF - Portable Document Format PLT - PLT / HPGL Format SHP - ESRI Shape Format SVG - Scalable Vector Graphics Format VEC - LEADTOOLS Vector Dump Format VPG - WordPerfect Vector Format WMF - Windows Metafile Formats WMZ - Compressed Windows Metafile Format </p>	
Video	<p> AVI DVD (.ifo, .vob) (requires a DVD-compatible video decoder Also known as a decompressor Also known as a decoder, this is a module or algorithm to decompress data. , this is a module or algorithm to decompress data.decompressor , this is a module or algorithm to decompress data. . For more information, refer to VOB Files .</p> <p> Mpeg1 Layer &ndash; 3 Audio (MP3) MPEG1 MPEG2 (requires an MPEG2-compatible video decoder such as LEAD MPEG2 Decoder) Ogg Quicktime WAV Windows Media Formats (WMA, WMV, ASF) DICOM ISO (MP4, AAC, AVC) Flash (FLV, F4V) </p>	Stöds ej.