



Big and easy to understand

GEWA Maxi

GEWA Maxi is a big, robust and easy-to-use remote control, which helps you to manage the equipment in your house. You can control your TV, sound system and lights, as well as open and close windows and doors.

If you have impaired mobility, vision, or if you struggle to understand a normal remote control, GEWA Maxi will guide you with pictures and buttons.

GEWA Maxi allows you to be more independent. You can make your own decisions, and not rely on help from family members or other carers.

Art.no 461130



ABILIA

GEWA MAXI MAKES DIFFICULT TASKS EASIER

Traditional remote controls often have too many buttons that are too small. They can be difficult to see, understand and press. Sometimes, you need more than one remote to watch TV, e.g. if you have a cable-box connected to the TV.

GEWA Maxi gives you a limited number of big buttons, that only show what you are wanting to control. Moreover, the press of a button can generate a series of events, e.g. first you switch on the cable-box and then the TV.



EASY TO CONTROL

GEWA Maxi has 16 big buttons that only require a little bit of pressure.

The buttons are lowered to facilitate that you hit the right button if you have tactile or motor-skill impairment.

Confirmation in the form of sound, light and vibration provides feedback that you have pressed a button.

If your hands shake, e.g. due to tremors, there is a possibility for activation delay. This means that a button must continue to be pressed for a selected amount of time before the choice is activated.

EASY TO UNDERSTAND

The buttons are adapted using supportive pictures that clearly symbolise what system the button controls.

Picture overlays are easily created via Abilia's website. A plastic sleeve protects the overlay against dirt and fluids.

If you do not need to use all the 16 buttons, you can hide the unused buttons using plastic caps. The plastic caps makes the surface feel smooth, and the button will not be visible.

EASY TO MANAGE

You can start using your GEWA Maxi quickly. A guide will show you how to e.g. transfer channels from the TV's usual remote control.

GEWA Maxi can be fastened to a table by using the included velcro strip. This will prevent it from sliding off or getting lost. In contrast to a traditional remote control, you do not need to aim GEWA Maxi at the equipment you wish to control.



TECHNICAL INFORMATION		ACCESSORIES		DELIVERS WITH
Length	242 mm	Product	Artikelnr.	Plastic sleeve
Width	160 mm	Plastic sleeve GEWA Maxi	461172	Velcro strips
Front height	23 mm	Velcro strips GEWA Maxi	461173	Plastic caps
Back height	42 mm	Plastic caps	461175	Battery
Weight	700 g			Manual