

A clear view of the day


MEMO Dayboard

MEMO Dayboard is a planning tool for those of you who need a visual overview of the day in relation to the course of time.

MEMO Dayboard is used to increase predictability in everyday life. Knowing what is going to happen, as well as when it will happen, is a basic human need.

MEMO Dayboard consists of a magnetic whiteboard and two time pillars, a red one for the day and a blue one for the night. The different colours make it easy to see whether it is day or night. By having text and pictures on the board, you can see what will happen during the day. You can also receive a signal that reminds you of when an activity should take place.

Art.nr: 508820 MEMO Dayboard, 12 hr
Art.nr: 508825 MEMO Dayboard, 24 hr



ABILIA


SEE WHAT WILL HAPPEN

Show the day's events and activities by writing or drawing directly on the board, or by adding pictures.



PRINT PICTURES FROM ABILIA PRINT-IT

Via the Abilia Print-It web service, you can easily print out Handi-pictures in a suitable size.



USE MAGNETS

Magnets make it easy to put up the information you want on the board. For example, the accessory contains Customisation Materials weekday magnets and holders for your own pictures.

HELP WITH TIME MANAGEMENT

The board's time pillars help you to keep track of the time of day. The dots on the time pillars go out as the day or night passes. By placing events and activities against the time pillar, you can see when they start or end.

BLUE LIGHT AT NIGHT

At night, the blue time pillar lights up and makes it easy for you to see that it is night. Studies show that this feature helps both you and your family members to sleep more calmly.



GET REMINDERS

If you need to be aware that something will happen at a certain time, you can set a reminder in the form of an audible signal. To make it easier to sign the reminder, you can connect an external switch.



SIMPLIFY WITH 12-HOUR DISPLAY

MEMO Dayboard is available in two models, where the time pillar either has a 24-hour or 12-hour display. The 12-hour model may be easier to use, as it corresponds to daily numbers.

TECHNICAL INFORMATION

Size	713 x 392 x 32 mm
Weight	3,8 kg
Power cord	9 V DC, 1A

PRODUCT

MEMO Dayboard - 12-hour version
MEMO Dayboard - 24-hour version

ART NR

508820
508825

ACCESSORIES

Customisation material	508785
Table stand	508794
Lockable door	508796
Pen Holder	508795

MEMO Dayboard is a CE-marked medical technology device.

