



Sense of time

MEMO Timer

MEMO Timer gives a visual representation of time, making something which otherwise is invisible, visible! Time is an abstract concept so learning it can be difficult for many people, when someone says half an hour, what does that mean? How do we feel that half an hour has passed? MEMO Timer allows you to get a quick and easy visual track of the time, allowing it to become accessible.

MEMO Timer displays the time as LED's that go out as the time passes, when the timer is finished the display LEDs will flash; this can be accompanied by an auditory signal or a tactile vibrate option.

MEMO Timer is waterproof so it can be used in the bathroom or outside, whatever the weather.

Art.no 500155



ABILIA

TAKE CONTROL OF YOUR DAY

Routines can be stressful and cause anxiety when you find understanding time difficult. MEMO Timer can be a helpful reminder in many ways, if you have to leave the house for school or work in 60 minutes, a 60 minute timer can be used to show how much time has passed and how much is left. The sense of control minimises the risk of unnecessary stress and anxiety.



The MEMO Timer can help you stay on track by being a reminder to get out the shower or bath after a selected period of time, this again allows you to relax knowing a reminder is in place when time is up rather than relying on another person to let you know.

PORTABLE AND EASY TO USE



MEMO Timer is easy to use and understand. It gives a clear picture of the amount of time thanks to the clear LEDs that periodically disappear as time passes by, which gives a visual impression of how much time remains. With support from light, sound and vibration, you will be informed when the selected time has come to an end.

MEMO Timer is small and smooth enough to fit into a normal sized pocket and is comfortable to hold in your hand. It's robust and

MEMO Timer is available in four versions:

MEMO Timer 8 shows time in 2, 4, 6, 8 minutes.
MEMO Timer 20 shows time in 5, 10, 15, 20 minutes.
MEMO Timer 60 shows time in 15, 30, 45, 60 minutes.
MEMO Timer 80 shows time in 20, 40, 60, 80 minutes.

PRODUCT

MEMO Timer 8	500150
MEMO Timer 20	500155
MEMO Timer 60	500160
MEMO Timer 80	500165

ART.NO

TECHNICAL INFORMATION

Size:	107x46x25 mm
Weight:	0.12 kg
Battery:	AA batteries 1.5 V
Volume:	max 80 dB
IP code:	IP67

MEMO TIMER COMES WITH

- MEMO Timer
- Wrist strap
- Neck strap
- 2x AA batteries 1.5 V
- Manual