



## Helps to keep track of the day

MEMOday

MEMOday is an aid for people who need extra support in keeping track of the time of day.

MEMOday shows easy-to-read information about the date, time, and whether it is morning, daytime, evening or night. One can also have the information read out loud.

For people who struggle to keep track of the time of day, the time support aid, MEMOday, provides an increased feeling of control, and thus a reduced sense of stress and anxiety. By showing whether it is night or day, MEMOday can, for example, prevent the person from going out in the middle of the night, or phoning family members at unsuitable times.

Art.nr: 508270



**ABILIA**

## EASY TO USE

With its stylish and attractive design, MEMOday blends right in with the decor of your home.

Place MEMOday in an area where you spend a lot of time. For example, you could place it in the bedroom, at your TV, or in the kitchen.

The information shown on the screen is easy to read. The text is large and has a high contrast. You can adjust the screen to show black text against a white background, or yellow text against a black background. At night, the screen automatically becomes darker, in order to not disturb your sleep.

Many people are helped by receiving information via their sense of hearing. Therefore, MEMOday has speech support that reads out the information when you press on the screen.



## EASY TO ADJUST

It is easy to set up MEMOday, and you can adjust the information that is shown on the screen according to your own wishes.

You can select an analogue or digital clock, or both of them at the same time. You can also raise and lower both the sound and the talking speed of the speech support.

You can reduce the amount of information on the screen, in order to make MEMOday even simpler and clearer.



In the event of a power failure, the text will disappear, but MEMOday will function as usual once the power comes back on.

For many people, it is important to receive information in their mother-tongue. MEMOday is currently available in 28 different languages.

MEMOday is a CE marked medical device.

### TECHNICAL SPECIFICATION

Dimensions (HxWxD)	140 x 185 x 95 mm.
Weight	0.35 kg.
Display	7 inch.
Available in following languages	English, Swedish, Norwegian, Danish, Finnish, Sami, Icelandic*, German, Dutch, Flemish, French, Spanish, Italian, Slovenian*, Bosnian*, Croatian*, Serbian*, Czech, Polish, Albanian*, Arabic, Persian*, Turkish, Somali*, Hungarian*, Russian, Greek, Cannada. (*= display reader not available)