

Art.nr: 464000

Emfit epilepsy alarm is designed for those who have epilepsy.

The alarm registers tonic-clonic seizures while you're lying in your bed.





EPILEPSY ALARM

The Emfit Epilepsy Alarm consists of a sensor that is placed beneath the mattress. The sensor is connected to a control device which can be mounted onto a wall using the accompanying bracket, or placed on a bedside table.

The sensor sends a signal to the control unit which sounds an adjustable alarm. It is also possible for the control unit to send alarm to a further alarm, such as pager or cellphone using our accessories All in One.

The signal system gives you the security you need as it calls for attention from family members or staff. The technical solution is adapted to the property's physical environment and the procedures for how the alarm is received.

Emfit epilepsy alarms are part of a series of reliable alarm products which are well tested and used in both home and residential care around the world.



TECHNICAL INFORMATION

SENSOR

Size 400 x 580 x 1,5 mm

Material Polyester

Cabel lenght 300 mm

CONTROL DEVICE

Size 97 x 127 x 30 mm Weight 150 g Color White

Battery 2 st 1,5 V, AA-battery

ART.NO ACCESSORIES

464015 Power adapter 464005 Wall mount 330701 All in One* 330700 All in One Basic*

* Not available at all markets

