

Small but strong IR-15 Wrist Transmitter

IR-1S Wristband transmitter is an IR remote control with only one button intended to assist persons who are unable to use a standard remote control due to motor, tactile or visual impairments.

The device can help the user to control electronic equipment at home, such as lighting and door openers. The device is intended for indoor and outdoor use. The bracelet transmitter can be used with the included bracelet or necklace.

IR-1S Wristband transmitter is a 1-channel IR transmitter that can be programmed with 7 different channels.

The product is marketed as a technical aid class I for people with disabilities and complies with Regulation (EU) 2017/745 - MDR.

Accessories:

Technical information:

Physical goals: (L x B x H): 46 x 44 x 14 mm Battery pack, art nr 401930 Weight incl batteries: 25 gram





Art nr 400500