



Feel more confident Emfit SafeFloor

SafeFloor is specifically designed for you who should not leave your bed or apartment unassisted.

With SafeFloor you will feel more secure at home. The system reduces the risk of you wandering around in the middle of the night. The system will alarm your nurse or relatives if you fall out of, or leave the bed.

Art.nr: 464020



ABILIA

SAFEFLOOR

SafeFloor contains one control device and one floor sensor.

When used at a bed-side, the SafeFloor enables an early warning for fall prevention as it gives the notification of the moment your feet touch the sensor. This can save many critical seconds, allowing the nurse to better prevent accidental falls.

CONNECT THE SYSTEM

The SafeFloor monitor has an audible alarm with adjustable volume and dry-contact output for connection to most nurse call systems or personal emergency phones.

The sensor sends a signal to the control unit which sounds an adjustable alarm that can be used to alert nearby nursing staff. It is also possible for the control unit to send to a further alarm, such as pager or cell-phone, using our accessories All in One.

Emfit Safe Floor are part of a series of reliable alarm products and are used in both home and residential care around the world.



SafeFloor enables an early warning the moment your feet touch the floor.

TECHNICAL INFORMATION

SENSOR

Size 927 x 580 x 2,4 mm
Weight 1,6 kg
Color Grey
Cable 300 mm

CONTROL DEVICE

Size 97 x 127 x 30 mm
Weight 150 g
Color White
Battery 2pcs 1.5 V AA battery

ART.NR ACCESSORIES

464015 Power supply
464024 5m extension cable
464005 Wall mount
330701 All in One*
330700 All in One Basic*

* Contact your Abilia representative to see if SafeFloor is available in your country.

