



Manage many functions

Control 18 & Prog - Remote controls

Control 18 and Prog has a large number of buttons that you can program, making it possible for you to manage various functions within your home.

The remote controls are intuitive and easy to use, and are suitable for the elderly or people who suffer from any form of cognitive or motor skills disability, or for those who cannot use an ordinary remote control for other reasons. Our remote controls provide you with a complete solution, which allows you to manage everything using the same remote control, while allowing you to clearly see which function you are controlling.

The ability to control functions in your surroundings yourself can give you increased independence and integrity in your home.



MANAGE ALL FUNCTIONS USING THE SAME REMOTE CONTROL

By using Control remote controls, you can manage all the functions in your home using the same remote control. You no longer need to learn or keep track of several different remote controls. You can use the remote control outdoors, to open entrance doors, garage doors, etc. When indoors, you can use the remote control for e.g. elevators, your front door, the lighting, radio, TV, stereo and the alarm.

The remote controls are programmable and can learn IR codes from virtually all remote controls that are available on the market, regardless of manufacturer or model.



Controller Overlay Keypad

EASY TO USE

The remote controls have large buttons that require very little pressure, which make it easy for you to push the correct button. When pressed, they provide feedback via audio and light signal.

You choose the functions in your home that you want to manage, and you label the buttons with the symbols, pictures or text that you need to remember and understand the various functions.

You can create your own labelling system, either by creating an overlay using the Overlay Design Web tool, or by using Word templates. Both of these are available on Abilia's website and contain completed symbols. You place your overlay over the buttons of your remote control.

CONTROL A SERIES OF EVENTS AUTOMATICALLY

You can create sequences that allow a single button to replace an entire series of button pushes and automatically perform a series of events. This makes it easier to e.g. find your favourite programme or a film channel. An additional application could be that the sound on the TV switches off when you answer your door's intercom.

The remote controls are able to transmit a uniquely coded channel for e.g. a security gate or front door. This provides you with increased security.

CHOOSE BETWEEN TWO DIFFERENT MODELS

Control 18 - More buttons provide more function possibilities

The remote control has 18 buttons with the possibility of performing 81 different functions.



Control Prog – Have control with your own control buttons

The remote control has 18 buttons with the possibility of performing 241 different functions.



Beyond the functions mentioned above, Control Prog has built-in scanning, which is convenient if you struggle to press the buttons directly. You are able to use external switches to manage all the functions.

COMES WITH

A changeable overlay with numbers, fingerguide, batteries LR03, neck strap, and manual.

TECHNICAL INFORMATION		PRODUCT	ART NR	ACCESSORIES	ART NR
Size	134 x 70 x 28 mm				
Weight	180 g	Control 18	425718	Control 18/Prog - Sheet of button labels	425776
Pressure	160 g	Control Prog	425700	Wall mount	810510
Pressure Control Prog	100 g				