



To keep track of the time of day

MEMOday

MEMOday is a tool for people who need extra support to keep track of the time of day.

MEMOday displays easy-to-read information about the date, time and whether it is morning, afternoon, evening or night. The information can also be read aloud.

For people who find it difficult to keep track of the time of day, the time support feature in MEMOday provides a greater sense of control and thus reduces stress and anxiety. By showing whether it is night or day, MEMOday can, for example, prevent the person from going out in the middle of the night or calling relatives at inappropriate times.

Art.nr: 508290



ABILIA


How to use MEMOday

With its stylish and attractive design, MEMOday blends seamlessly into the home.

Place MEMOday in a spot where the user spends a lot of time. This could be in the bedroom, by the TV or in the kitchen, for example.

The information displayed on the screen is easy to read. The text is large and high-contrast. At night, the screen automatically dims so as not to disturb sleep.

Many people find it helpful to take in information through their sense of hearing. That is why MEMOday has a text-to-speech feature that reads out the information when the user taps the screen.



Customising MEMOday

It is easy to set up MEMOday with step-by-step instructions. The information displayed on the screen can be customised to suit the user's needs and preferences.

You can display an analogue or digital clock, or both at the same time. You can adjust both the volume and the speech rate of the speech support.

The amount of information on the screen can be reduced to make MEMOday even simpler and clearer.

More great features

In the event of a power cut, the text disappears, but MEMOday works as normal once the power is restored.

For many people, it is important to receive information in their mother tongue. MEMOday is currently available in 28 different languages, which can be easily selected in the settings menu.

MEMOday is a Class I medical device for people with disabilities and complies with Regulation (EU) 2017/745 – the MDR.

TECHNICAL SPECIFICATION

Dimensions (HxWxD)	140 x 185 x 95 mm
Weight	350 g
Screen	7 inches
Available languages	Swedish, Norwegian, Danish, Finnish, Icelandic*, Sami, English, German, Dutch, Flemish, French, Spanish, Italian, Slovenian*, Bosnian*, Croatian*, Serbian*, Czech, Polish, Albanian*, Arabic, Persian*, Turkish, Somali*, Hungarian*, Russian, Greek, Kannada*. (Languages marked with * do not have voice support)