

## Gewa Control Medi

- A tiny IR transmitter with advanced features



Gewa Control Medi


*Facing slight of hand*

*Easy to press buttons*

*Can be customized with image support*

*Ability to use scanning*





## CONTROL MEDI

Control Medi is a newly developed IR transmitter with new design. Small and compact yet with large, easy to press buttons, Control Medi is created for people with motor disabilities.

Control Medi is available in three different models. The most advanced feature has 8 buttons and the ability to send radio signals (GewaRadio). It also has a built-in scanning option for people who cannot press the buttons directly. The simplest model has 4 or 2 buttons, two keyguards are included.

Naturally, all the models can learn IR-signals from virtually all remote controls in the world.

When you press a button, confirmation is given as vibration, sound and light.

It is possible to go to Abilia's website and create your own button labels.

## MODELS AND FEATURES

### 429200 Control Medi Multi

- 8 buttons
- Scanning
- GewaRadio
- IR learning

### 429220 Control Medi Standard

- 8 buttons
- IR learning

### 429240 Control Medi Easy

- 4 or 2 buttons
- IR learning

## SPECIFICATIONS

Size of casing:	56 x 100 x 22 mm
Size (8) buttons:	15 x 15 mm
Size (4/2) buttons:	15 x 36 mm
Pressure force:	< 160 g

## MODELS

> **ITEM:** 429200 Control Medi Multi

---

> **ITEM:** 429220 Control Medi Standard

---

> **ITEM:** 429240 Control Medi Easy

---

> **MORE INFORMATION:** Visit our website for more photos, videos and information

---