


ALARM

All in One

A universal solution for forwarding alarms and alerting colleagues



All in One sends an alert that help is needed. Registers alarms and forwards them to the receiver or a mobile phone.



ALL IN ONE

All in One is a universal solution for forwarding alarms. The solution consists of two main units – a transmitter and a receiver, often referred to as a pager.

All in One can forward alarms from all types of wired sensors and a number of wireless transmitters, such as transmitter bracelets and universal transmitters.

All in One can forward alarms to the receiver as well as to a mobile phone.

Typical applications include forwarding triggered epilepsy alarms or alerts that a colleague needs help.

SAFETY

Great focus was placed on ensuring safety and a good range. All in One has a so-called heartbeat function that alerts you if you go out of range, or if there is no contact between the transmitter and the receiver for some reason.

ALL IN ONE

- Good range, over 600 m if unobstructed
- Possible to send an alarm to both the receiver and a mobile phone
- Possible to forward alarms from both wired and wireless sensors
- Heartbeat function, safety function that alerts you if there is no contact between the transmitter and receiver
- Intuitive, simple installation and easy to use
- The receiver can receive alarms from 4 different universal transmitters

ALL IN ONE BASIC

All in One is also available in a simpler version – All in One Basic. All in One Basic can forward alarms from wired sensors to a hand-held receiver.

All in One Basic has the same good range and safety as the main model.

All in One

> ART NR: 330701 All in One

> ART NR: 330700 All in One Basic

> ART NR: 330706 Motion sensor

> ART NR: 330708 Universal transmitter

> ART NR: 330709 Transmitter bracelet

> ART NR: 330720 External Vibrator

> VISIT OUR WEBSITE for more information