



For a safe rest

Emfit SafeBed

Art.nr: 464001

Emfit SafeBed is a "missing person alarm" that can register when a person gets out of bed, or fails to return to bed within a pre-set time.

This provides safety for you if you have a tendency to rise up and walk away unattended as the alarm draws attention of personnel or relatives who are in the vicinity.



ABILIA

SAFE BED

Emfit SafeBed consists of a sensor that is placed under the mattress, and a main unit where settings are entered.

If the time setting is short (ie. 3 seconds) the alarm will sound when the person gets out of bed. By increasing the time (up to 30 minutes) the alarm will sound if the person does not return to bed within the pre-set time.

The sensor sends a signal to the control unit which sounds an adjustable alarm that can be used to alert nearby nursing staff. It is also possible for the control unit to send to a further alarm, such as pager or cellphone, using our accessories All in One.

Emfit epilepsy alarms are part of a series of reliable alarm products which are well tested and used in both home and residential care around the world.



The sensor is placed under the mattress

TECHNICAL INFORMATION

SENSOR

Size 400 x 580 x 1,5 mm
 Material Polyester
 Cable length 300 mm

CONTROL DEVICE

Size 97 x 127 x 30 mm
 Weight 150 g
 Color White
 Battery 2pcs 1,5 V, AA-battery

ART.NO ACCESSORIES

464015 Power supply
 464005 Wall mount
 330701 All in One*
 330700 All in One Basic*

* Not available at all markets