



For a social life
Lightwriter® series

The Lightwriter SL40 is a dedicated device for those who are literate, but unable to communicate effectively through speech. With a Lightwriter you can participate in social interactions while at work or leisure, even in noisy, sunny outdoor or poorly lit environments.

Using dual screens, your communication partner can read the text whilst maintaining a natural face-to-face dialogue, including engagement through gestures and body language.



ABILIA




The Lightwriter SL40 is a dedicated device for those who are literate, but unable to communicate effectively through speech. With a Lightwriter you can participate in social interactions while at work or leisure, even in noisy, sunny outdoor or poorly lit environments.

FOR A MORE ACTIVE LIFE

Do you want to socialize with your friends and communicate with them in different environments and settings? The Lightwriter gives you access to technology specifically designed for communication, that can help enrich your life and encourage and support engagement with your family, friends and the community.

STAY IN TOUCH WITH YOUR FRIENDS AND FAMILY

To stay in touch with friends and family is important for everyone. The Lightwriter lets you send text messages with a SIM card inserted. With the SL40 Connect you can also make phone calls and back-up your personalised vocabulary to a USB/ Personal computer. In the notebook, you can save important textmessages and your own notes.



EASY COMMUNICATION

Lightwriter has two bright screens, so noise and poor lighting is no obstacle for clear communication. Using dual screens, your communication partner can read the text whilst maintaining a natural face-to-face dialogue, including engagement through gestures and body language.

The displays are designed to be easily readable, even outdoors in the sunlight. The built-in speaker has very good sound quality and faces your communication partner. The Lightwriter is easy to use, press a button and you are all set!



"I like that I can talk to many different people when I'm out"

BRING IT EVERYWHERE

The Lightwriter is small, portable and robust. It comes with a carrying case. You can take it everywhere as it is designed to withstand the bumps and knocks of everyday life. It is possible to add a protective cover which shields your Lightwriter keys from splashes of liquid (available separately). Long battery life means you can use your Lightwriter as much as you want all day every day, if you recharge it every night.

EASY TO USE

Lightwriter has a user friendly menu and excellent sound quality using text-to-speech technology (Acapela). It can predict which words you are going to use and proposes various alternatives. In addition to writing text it is possible to pre-store words and sentences for even faster communication.

Each Lightwriter comes with a set of keyguards. Keyguards can help to guide your fingers to the correct location and avoid accidental key presses, ideal if you have poor fine motor control.

OTHER FEATURES THAT SIMPLIFY YOUR LIFE

There are several additional features included in the Lightwriter that will simplify your life, such as alarms to remind you of important events. You can also use the Lightwriter as an environmental control and as a remote control for your TV.



WANT TO KNOW MORE?

You will find additional features, other than those described in this leaflet, in the user guide that can be found on our website.

ABOUT ABILIA

Abilia develops, manufactures and provides assistive technology for those who need support in the areas of Communication, Cognition and Environment Control & Alarm.

CONTACT US

Phone +46 (0)8-594 694 00

E-mail: info@abilia.se

Visit www.abilia.com for more information about our assistive devices and solutions.



PRODUCT MODELS

Lightwriter SL40
Lightwriter SL40 Connect

LIGHTWRITER COMES COMPLETE WITH

Three different keyguards, carry bag, charger, tool kit, user guide and getting started leaflet, memory stick (Connect model only).

ACCESSORIES

In Use Bag
Trabasack Mini - Slim Lap Tray Bag
Complete Hand, Neck & Waist Strap Set
In-car charger

TECHNICAL SPECIFICATIONS

Size:	240mm x 140mm x 55mm
Weight:	910 g
Speech synthesis:	Acapela HQ
Battery:	Rechargeable nickel metal hydride
Battery life:	All day battery life with standard usage