



## Speak with confidence

### Micro Voice Amp

The Micro Voice Amp is perfect for those with a weak or strained voice, or who regularly use their voice intensively.

Whether you are healing from surgery or have more permanent damage to your vocal chords, the Micro Voice Amp is designed to allow you to get your message across loud and clear without straining or causing further damage to your voice.

Micro Voice Amp gives you control over the volume of your voice in noisy environments, to groups of people, indoors and outdoors with confidence.



**ABILIA**  
Toby Churchill

## ALL YOU WANT FROM A VOICE AMP

Based on your feedback we have improved the sound quality, giving you a loud, clear voice with reduced distortion and feedback sounds.

The battery is not only fast at charging, but now has an extremely long life so you can forget about worrying if it has been charged every time you go out.

## EASY TO USE

We have only included the features you need, and kept operation intuitive and simple - just switch it on and start talking into the microphone. The coloured LED indicators help you to understand what mode you are in.



## DESIGNED FOR DAILY LIFE


Designed with a discreet compact black case, the Micro Voice Amp is lightweight and portable enough to wear or carry.

The Mute function allows you to continue to talk privately to those close to you without having to turn off your amplifier.



## CUSTOMISE TO SUIT YOUR LIFESTYLE

Choose from a high quality AVI headset microphone; a discreet lapel microphone which can be clipped onto clothing; or a transdermal microphone which is fitted near your throat - designed to pick up even faint vocal sounds.



Lapel Clip



AVI Headset



Transdermal

It is also possible to use the audio output to plug into an external device such as a loudspeaker.

MICRO VOICE AMP COMES COMPLETE WITH  
Charger, belt clip and waist strap, user guide,  
one microphone option.

## TECHNICAL SPECIFICATIONS

Size: 122 x 67 x 32mm  
Approx Weight: 210g  
Ports: 3.5mm mono input,  
3.5mm stereo output.  
Battery: Integrated rechargeable battery  
40hrs active, 200hrs standby.

## MICROPHONE OPTIONS

