

# RELEASE NOTES - ROLLTALK DESIGNER 8.4.2

## **Memoplanner support**

Memoplanner is now running on the current version of Rolltalk Designer.

## **2965: Nova: Volume not set when changing TTS engine**

Update default volume when changing voice technology

## **3015: Album slideshow does not work**

Updated Rolltalk Designer - Examples Page Library Camera and Field Library Album

## **3165: Runtime language inadvertently set equal to master language**

Fixed program internal issue.

## **3136: Crash when connected to Multibox on Win10 AU**

Fixed problem with connection to Multibox on Windows 10 Anniversary Update.